

HTM2010 Porous Load - Yearly Tests

Porous Load Yearly Summary					
Rev. 1.3.0011					
Authority					
Hospital					
Department	Area R5				
Date of Report	08/01/2010 12:21:55	Dates of Tests		05/01/2010	
Steriliser	Manufacturer		L.T.E. Scientific Ltd		
	Serial Number		L7096		
	Reference number		N/A		
	Model		Touchclave K150E Porous Load		
	Usable Chamber Space		150Litre	litres	
	Date of Manufacture		18/09/06		
Service Provider	L.T.E. Scientific Ltd				
Results of Yearly Tests					
Tests as Specified in HT2010	Ref	Pass/Fail	Cycle No.	Result	
Safety checks	5.8	Pass	N/A		
Steam non-condensable gas test	9.4	Pass	1933-1935	Concentration of NCG	1.09 %
Steam superheat test	9.20	Pass	1936	Superheat	0.23 °C
Steam dryness test	9.30	Pass	1937	Dryness Value	1.00
Vacuum leak test (before sensors)	11.2	Pass	1938	Leak rate	0.2 mbar/min
Vacuum leak test (sensors inserted)	11.2	Pass	1940	Leak rate	1.0 mbar/min
Automatic control test	12.1	Pass	1941+1942	ST selected 134&12°C	ST hold time 4&20Mins
Verification of calibration	12.2	Pass	1941+1942	See below	
Air detector performance test small load	11.45	Pass	1945	Leak rate	9.7 mbar/min
Air detector performance test full load	11.53	Pass	1946	Leak rate	9.7 mbar/min
Thermometric test full load	13.15	Pass	1949+1952	See below	
Thermometric test small load	13.7	Pass	1948+1951	See below	
Performance Re-Qualification	8.64	N/A	N/A	See Performance Re-Qualification Test (if applicable).	
Vacuum leak test (sensors removed)	11.2	Pass	1953	Leak rate	0.2 mbar/min
Air detector function test	11.60	Pass	1955	Air Detector Setting	9.7 mbar
Bowie-Dick test for steam penetration	13.39	Pass	1956	Type of Test Pack	Browne TST

List of Calibrated Equipment					
Data Logger			Heat Source		
Manufacturer	Fluke		Manufacturer	Isotech	
Model	NetDAQ		Model	Fast-Cal	
Serial No.	LTE-LOG16MB		Serial No.	LT343	
Calibration House	N/A		Calibration House	TMS Brockwell	
Cert No.	N/A		Cert No.	13802	
Renew date	N/A		Renew date	05/10/10	
Temp Ref Probe			Pressure Calibrator		
Manufacturer	Isotech		Manufacturer	GE Druck	
Model	Fast-Cal Indicator		Model	DP1103	
Serial No.	LT343		Serial No.	DPVTG5	
Calibration House	TMS Brockwell		Calibration House	TMS Brockwell Calibration	
Cert No.	13802		Cert No.	12773	
Renew date	05/10/10		Renew date	03/03/10	

HTM2010 Porous Load - Yearly Tests

TQReports
WWW.TQREPORTS.INFO

Yearly Safety Checks		
Safety Check	Completed	Comments
Safety valve(s) operation	Yes	
Condition of door seal	Yes	Replaced Door Seal.
Door interlocking mechanism	Yes	
Door interlock pressure operation	Yes	
Door steam generation / admission interlock	Yes	
Housekeeping	Yes	
Daily record check	N/A	
Weekly record check	N/A	
Quarterly record check	N/A	
Yearly maintenance	Yes	
	N/A	
	N/A	
Results	Pass	

HTM 2010 - 5.8

Steam Non-Condensable Gas Test Porous Load						
	Test 1		Test 2		Test 3	
Cycle No.	1933		1934		1935	
Volume of Gas Collected (Vb)	0.8 cubic cm		0.9 cubic cm		1.8 cubic cm	
Vol. of Condensate Collected (Vc)	105 cubic cm		110 cubic cm		105 cubic cm	
Fraction of Non-Condensable Gases: 100 x (Vb / Vc) =	Yearly 0.76 %	Validation N/A %	Yearly 0.82 %	Validation N/A %	Yearly 1.71 %	Validation N/A %
Result	Pass		Pass		Pass	
Comments	Condensable gases should not exceed 3.5%.		Condensable gases should not exceed 3.5%.		Condensable gases should not exceed 3.5%.	

HTM 2010 - 9.4

Steam Superheat Test			
Cycle No.	1936		
Temperature in expansion Tube (Te)	99.69 °C	Superheat = Te - To =	0.23 °C
Temp. of boiling Water at local atmos. Pressure (To)	99.46 °C	Validation	N/A
Steam service supply Temperature	110.49 °C	Result	Pass

HTM 2010 - 9.20

Steam Dryness Test			
Cycle No.	1937		
Weight of flask assembly when empty		M1 =	816.5 g
Weight of flask assembly + cold water		M2 =	1443.3 g
Initial weight of water in flask		Mw = M2 - M1 =	626.80 g
Initial temperature of water in flask		To =	21.86 °C
Average temperature of steam delivered to Sterilizer		Ts =	139.43 °C
Final temperature of water and condensate in flask		T1 =	80.71 °C
Weight of flask assembly + condensate collected		M3 =	1513.5 g
Weight of condensate collected		Mc = M3 - M2 =	70.20 g
Latent heat of dry saturated steam at temperature Ts		L =	2145.54 kJ/kg
Dryness value (T1 - To) (4.18Mw + 0.24) / (LMc) - 4.18(Ts - T1) / L		D =	1.00 Validation N/A
Result			Pass

HTM 2010 - 9.30

HTM2010 Porous Load - Yearly Tests

WWW.TOPREPORTS.INFO

Rev. 1.2.0010 Porous Load Automatic Control Test					
Cycle No.	1941	Sterilisation Temp	134.0 °C	Yearly	Validation
Air Removal	Negative Pulsing	Minimum Vacuum		306 mBA	N/A
		Maximum Vacuum		105 mBA	N/A
		Number		5	N/A
		Duration		00:01:45 hh:mm:ss	N/A hh:mm:ss
	Positive Pulsing	Minimum Pressure		1192 mBA	N/A
		Maximum Pressure		2020 mBA	N/A
		Number		5	N/A
		Duration		00:02:07 hh:mm:ss	N/A hh:mm:ss
Steam Admission and Sterilizing	Time to Attain Sterilising		00:06:55 hh:mm:ss	N/A hh:mm:ss	
	Heat Up Time		00:01:46 hh:mm:ss	N/A hh:mm:ss	
	Indicated Pressure	Start		3200 mBA	N/A
		Mid		3244 mBA	N/A
		Max		3246 mBA	N/A
		End		3208 mBA	N/A
	Indicated Temperature	Start		135.8 °C	N/A
		Mid		136.1 °C	N/A
		Max		136.4 °C	N/A
		End		136.0 °C	N/A
	Recorded Pressure	Start		N/A mBA	N/A
		Mid		N/A mBA	N/A
		Max		N/A mBA	N/A
		End		N/A mBA	N/A
	Recorded Temperature	Start		135.3 °C	N/A
		Mid		135.9 °C	N/A
		Max		135.9 °C	N/A
		End		135.8 °C	N/A
	Measured Pressure	Start		3206 mBA	N/A
		Mid		3206 mBA	N/A
		Max		3245 mBA	N/A
		End		3218 mBA	N/A
	Measured Temperature	Start		136.96 °C	N/A
		Mid		136.38 °C	N/A
		Max		136.96 °C	N/A
		End		136.64 °C	N/A
	Steriliser Hold Time		00:04:00 hh:mm:ss	N/A hh:mm:ss	
	Drying	Maximum Vacuum		30 mBA	N/A
Duration		00:16:13 hh:mm:ss	N/A hh:mm:ss		
Air Admission	Duration		00:00:34 hh:mm:ss	N/A hh:mm:ss	
Total Cycle Time	Duration		00:27:42 hh:mm:ss	N/A hh:mm:ss	
				Result	Pass

HTM2010 - 12.1

HTM2010 Porous Load - Yearly Tests

Rev. 12.0010 Air Detector - Performance Test Small Load		
Cycle Number	Yearly	Validation
1945		
Load Temp. (Tp) at the Start of Plateau Period	133.32 °C	N/A
Drain Temp. (Tc) at the Start of Plateau Period	134.21 °C	N/A
Temperature Depression (Tc - Tp)	0.89 °C	N/A
Induced Leak Value	9.7 mbar/min	N/A
Air Detector Trigger Point Set Up	75.0 °C	N/A
Induced Leak to be used on future Air Detector Function Tests	9.7 mbar/min	N/A

HTM2010 - 11.45

Rev. 12.0010 Air Detector - Performance Test Full Load		
Cycle Number	Yearly	Validation
1946		
Load Temp. (Tp) at the Start of Plateau Period	132.76 °C	N/A
Drain Temp. (Tc) at the Start of Plateau Period	134.11 °C	N/A
Temperature Depression (Tc - Tp)	1.35 °C	N/A
Induced Leak Value	9.7 mbar/min	N/A
Air Detector Trigger Point Set Up	75.0 °C	N/A
Induced Leak to be used on future Air Detector Function Tests	9.7 mbar/min	N/A

HTM2010 - 11.53

TOReports
 WWW.TOREPORTS.INFO

HTM2010 Porous Load - Yearly Tests

TOReports
 WWW.TOREPORTS.INFO

Rev. 1.3.0013				Thermometric Test for Full Load	
Cycle No.	1949	Sterilisation Temp	134.0 °C	Yearly	Validation
Start of Plateau Period	Drain			134.02 °C	N/A
	Load			133.32 °C	N/A
	Top Sheet			133.44 °C	N/A
	Pressure			2953 mBA	N/A
Start of Holding Time	Drain			134.71 °C	N/A
	Load			134.08 °C	N/A
	Top Sheet			134.20 °C	N/A
	Pressure			3013 mBA	N/A
Mid of Holding Time	Drain			136.38 °C	N/A
	Load			136.30 °C	N/A
	Top Sheet			136.36 °C	N/A
	Pressure			3200 mBA	N/A
Max Values During Holding Time	Drain			136.79 °C	N/A
	Load			136.49 °C	N/A
	Top Sheet			136.56 °C	N/A
	Pressure			3217 mBA	N/A
End of Holding Time	Drain			134.24 °C	N/A
	Load			134.11 °C	N/A
	Top Sheet			134.07 °C	N/A
	Pressure			2983 mBA	N/A
Equilibration time				0:04 mm:ss	N/A mm:ss
Drain Fluctuation				0.21 °C	N/A
Load Fluctuation				0.26 °C	N/A
Max difference between Drain and Load				0.21 °C	N/A
Post vacuum stage time				00:19:51 hh:mm:ss	N/A hh:mm:ss
Pressure at end of vacuum hold				42 mBA	N/A
Load sensibly dry				Pass	N/A
Automatic test requirements met				Pass	N/A
Hold Time				00:04:19 hh:mm:ss	N/A hh:mm:ss
Total Cycle Time				00:41:41 hh:mm:ss	N/A hh:mm:ss
Result					Pass

HTM2010 - 13.15

Job Ref: Porous Load 134 4M
 Cycle No.: 1949
 Serial No.: L7096

Service Provider:
 Site Name:
 Operator:

L.T.E. Scientific Ltd

Test Name: HTM2010 Porous Load 134C
 Test Date: 07/01/2010
 Start Time: 11:31:37

