



# Dental

Equipment, Service & Training

LTE Scientific has over 60 years experience in the field of decontamination and sterilization. Our Dental division can deliver a comprehensive package of services, equipment, advice and training designed to meet the increasing demands of today's dental practices.

## LTE Dental

### Equipment

The LTE Mediclave range of dental autoclaves comes in 18 and 23-litre capacities and is available in B, S and N classes meaning that we have a model to suit all applications.

The separate steam generator keeps the chamber free from heaters, standing water and debris, thus reducing the potential risks of cross-contamination. With up to 13-pre-set cycles, Mediclave is one of the easiest machines to use and offers total flexibility. Contact our sales office for a brochure or download one today at [www.lte-scientific.co.uk](http://www.lte-scientific.co.uk).

In addition, LTE can provide a full package of equipment for the decontamination room. This includes a range of benchtop, under-bench and free-standing washer-disinfectors, ultrasonic baths and furniture specifically designed for the dental market.



18 Litre Model with Data Logger

## Service

LTE Service Centre employs a team of nationwide team of skilled service engineers, who are supported by technical engineers based at our head office.

Our engineers are certified to HTM 2010 & HTM 2030 and will be trained to HTM01 05 and HTM01 01 as these standards are introduced. LTE is also a UKAS accredited company.

In addition to supporting our own products, we test and validate all makes and types of decontamination equipment using the latest technology which enables the engineer to leave test results with the customer.

We ensure that large stocks of spare parts are always available and pride ourselves on the speed and quality of response that we are able to offer.

## HTM 01-05

Decontamination processes and procedures are an extremely important part of any dental facility, and the introduction of HTM 01-05 will bring about many changes for the dental market. LTE can assist in the maintenance and validation of your decontamination equipment, by providing cost effective maintenance and validation packages tailored to your requirements.

LTE has recently been selected for the NHS Supply Chain Dental Decontamination Framework agreement for maintenance and validation, giving you the confidence that that we have the capability to carry out this work.

LTE Service Centre can offer the following:

- Quarterly/Annual service & testing to HTM 01 -05
- Routine maintenance service contracts
- Installation & commissioning
- Emergency breakdown service
- Supply of spare parts
- In house repair service
- Training
- Application support & helpline

## Training and Support

The introduction of HTM01-05 will change the way in which dental practices need to approach decontamination both in terms of equipment, procedures and staff training.

Dental practices need to demonstrate that adequate training has been carried out in the area of decontamination. This is required to satisfy Continuing Professional Development as required by the General Dental Councils life-long learning scheme.

To meet this growing demand, we can offer a range of evidence based in-house and on-site training courses to suit individual client needs.

LTE employs a Clinical Nurse Advisor (CNA) who is a registered general nurse on part 1 of the register and has many years experience in the fields of decontamination, infection control and general/oral surgery in the NHS.

Essential quality requirements stipulate that certain functions need to be employed by dental practices in order to achieve the effective cleaning and inspection of instruments.

Our CNA can visit site to carry out an evaluation of your current decontamination & training provision, based on CE and regulatory compliance. We will then provide a written report outlining our findings and recommendations.

