



AE/MED/60

LARGE STEAM MEDICLAVE 60LT



A cost-effective, comprehensive autoclave for the sterilization of instruments, specifically geared to the dental, chiropody, surgical, and veterinary markets.

RELIABILITY

Fully compliant with all current directives

SECURITY

Automatic closing with double safety lock

“B” CLASS AUTOCLAVE

VACUUM PUMP

Provides optimal air removal for deep load penetration

LOW CONSUMPTION

Minimal water consumption per cycle



STERILIZATION CHAMBER CAPACITY

60lt

AUTOMATIC WATER LOADING AND DISCHARGE

BUILT-IN PRINTER

BENCH TOP
Ideal for areas with space constraints

TRACEABILITY

SD card slot as standard

TECHNICAL DATA

Power supply:	230V ~ 50/60Hz Monophase 15A
Chamber volume:	60lt in stainless steel (ø350x606mm)
Load:	9Kg
Weight:	62Kg
Dimensions without cabinet:	615 x 730 x 750 mm (WxDxH)
Dimensions with cabinet :	615 x 730 x 1500 mm (WxDxH)

Produced in accordance with the current directives and standards.

OPTIONAL EXTRAS

Floor Standing Cabinet
RO Water System

LTE is a progressive company and reserves the right to alter the specification of its products without notice. E&OE MEDL1 05-21