



# Eco-Drying Cabinets

ENERGY-EFFICIENT DRYING WITHOUT COMPROMISE



## Tackling Sustainability with our new Drying Cabinet!

Through recent innovation, our new Eco-Drying Cabinets can significantly reduce power usage compared to standard glassware drying cabinets.

Complete with digital temperature controller and electronic 7-day timer, the eco-drying cabinet range is designed to offer market leading operational efficiencies over standard drying cabinets while remaining economically priced.

### Features Include:

- Energy saving of up to 65% when compared to standard glassware drying cabinets
- Insulated on 5 sides
- Low heat output
- Digital PID temperature control
- Programmable 7-day timer

## TECHNICAL INFORMATION

Cat. No	Capacity (Litres)	Air Circulation	Max Temp (°C)	Internal Dims HWD (mm)	External Dims HWD (mm)	Shelves/ Positions
TD0100SN	100	Natural Convection	85	440 x 775 x 300	680 x 900 x 390	2/20
TD0180SN	180	Natural Convection	85	555 x 995 x 325	775 x 1125 x 420	2/34

### TYPICAL ENERGY USAGE AT 75 (°C)

Cat. No	Energy Used in 1hr (kWh)	Energy Used in 8hr (kWh)	Energy Used in 24 hr (kWh)
TD0100SN	0.201	1.608	4.824
TD0180SN	0.312	2.496	7.488

### HEAT UP TIMES & ENERGY USAGE (TD0100SN)

Set Temp	Load	Energy Used (kWh)	Time taken/ mins
22-75°C	No	0.193	20

Tests carried out at an ambient temperature of 22°C