Precision Cooling. The LTE Way

RELIABLE PERFORMANCE. PRECISE TEMPERATURE MANAGEMENT.

From everyday laboratory storage to temperature-critical applications, the new LTE refrigerator and freezer range delivers dependable solutions for protecting valuable samples and materials.

Designed for practical laboratory use, our new range of refrigerators and freezers provide precise, reliable and secure temperature-controlled storage for medicines, reagents and other temperature-sensitive products.

Our refrigerator range is available in three capacities, from 320 to 1,220 litres, alongside a 1,332-litre solid-door freezer.



KEY FEATURES:

- Stainless Steel interior which is easy to clean while remaining durable, supporting long-lasting performance (excluding TRS0130WD model)
- Sleek white external finish
- Visual and Audible Temperature Alarm for out-of-range temperatures, open door, power failure or sensor errors
- Safety Thermostat
- Touchscreen controller, enhancing user accessibility for easy of use
- Designed for Climate Class 4 which can operate in conditions up to 32°C and 60% humidity
- Adjustable shelving provides greater flexibility, allowing the internal space to accommodate products of different sizes
- Interior LED lighting, providing users with greater visibility
- Uniform cold distribution with automatic defrosting and defrost water evaporation

TECHNICAL INFORMATION

Model	Capacity (litres)	Temp Range °C	Inner Dims HWD (mm)	Overall Dims HWD (mm)	Weight (KG)	Shelves
TRS0320WS	320	+2°C—+8°C	1516 x 382 x 560	2005 x 490 x 716	102	5
TRP0320WS	320	+2°C—+8°C	1516 x 382 x 560	2005 x 490 x 716	102	5
TRS1220WS	1220	+2°C—+8°C	1510 x 1358 x 856	2037 x 1358 x 856	227	10
TFS1332WS	1332	-10°C- -26°C	1510 x 1203 x 680	2038 x 1356 x 857	227	10
TRS0130WD	130	+2°C—+8°C	616 x 452 x 504	850 x 626 x 600	45	3

EFFICIENCY AND SUSTAINABILITY

- Quiet operation with acoustic emissions <50 dB
- Energy rating B or A, reducing electrical consumption for improved sustainability.

OPTIONS AND ACCESSORIES

Description
Voltage: 110/60 Hz and 220/60 Hz
Stainless Steel shelves and Stainless Steel exterior finish
Pedal door opening for greater convenience

* Options are factory fit and need to be ordered at the same time as the cabinet.

LTE is a progressive company and reserves the right to alter the specification of its products without notice. E&OE

Need a different configuration or capacity? Talk to our team about the wider range of options available.

T: +44 (0) 1457 876221 **E:** info@lte-scientific.co.uk